



CATALOGUE DECUBITUS THERAPY

CLINIC RANGE



Products for more
sustainability!

Scan me!



Here you can find the
current online version.



CONTENTS

1. ABOUT FUNKE MEDICAL GMBH 5 – 10

Company information	7
Environmental responsibility	8
Certifications	9

2. PRODUCT FEATURES 11 – 14

Clinic range	13
Features of care mattresses	14

3. PRODUKTE 15 – 66

GREENLINE®	17
CLINIC CARE	23
EVAQ®	29
BARIATRIC CARE	37
AIRLINE®	43
MEMOLINE®	49
CLINIC CARE	51
EVAQ®	53
Positioning aids	57
Positioning pillows	61
Heel protector boot	81
Turning and positioning aid	83

4. SERVICES 87– 95





ABOUT FUNKE MEDICAL GMBH



Company information



Funke Medical GmbH – since 1988

The Funke Medical GmbH is a German company, which manufactures high quality medical systems for seating, lying and positioning of patients. The production takes place at three locations with a production area of 17.800m² in Raesfeld, in North Rhine-Westphalia. The products are manufactured exclusively in Germany.

The Funke Medical GmbH is the market leader in the field of decubitus prophylaxis and therapy for homecare in Germany. It is an international operating company, which already has long-standing business relationships with over 1.400 customers from the medical trade and clinics within Germany and abroad.

You would like to know more about our company and our product range? No problem!

Whether by e-mail, telephone or in person- the employees of Funke Medical GmbH are of course available to you at any time.

Further information can be found on our homepage klinikmatratze.de!

We look forward to hearing from you soon!

Dr. Ramona König (CEO), Helmut Wessels, Christian Maruhn



Administration and production facility



Warehouse and Service Centre



Administration and production facility

Environmental responsibility

Environmental and climate protection are very important to us. For this reason, we attach great importance to resource-saving processes in order to keep our ecological footprint as small as possible. To become more sustainable, we have taken the following measures:



Photovoltaic

Since 2012, we generate 250.000 kWh of electricity from solar energy using 4 photovoltaic systems. The produced amount covers about 65% of our own demand. This means that approx. 130 tons less CO₂ will be emitted into the environment every year.



Climate-friendly mattresses

Our **GREENLINE**® mattresses are made of **ECOLAST**® foam. Compared to conventional decubitus therapy mattresses (*standard foam mattresses or alternating pressure mattresses*) they have approx. 22%* less CO₂ emission.



Sustainable packaging

Our seat cushion cartons/HEEL PROTECT® cartons are made of grass paper. This is a fast-growing, sustainable raw material. The grass used for production comes from previously unused green areas, so-called compensatory areas. They are also only mown once at the beginning and once at the end of the year, so that the main vegetation periods of the grasses are not affected.



Water saving:

During the production of grass paper, significantly less water is required than for virgin wood pulp or the processing of waste paper.



CO₂ saving:

Grass paper produces significantly lower CO₂ emissions than virgin wood pulp.



Resource-saving:

Minimising the consumption of important resources such as wood.



Biodegradable:

Grass paper consists of natural fibres, which are biodegradable. The flexo inks used are water-soluble, harmless and can be recycled with the cardboard.



Saving of plastic

The seat cushions and HEEL PROTECT® are also provided with a fresh seal instead of an additional film. In this way we save 6 tons of foil per year.

We want to preserve our natural resources with all these measures in order to live up to our environmental responsibility.

THE WORLD IS A FINE PLACE
AND WORTH FIGHTING FOR.

Certifications

With the certification according to DIN EN ISO 13485 the Funke Medical GmbH is authorized to manufacture medical products and to sell them in accordance with relevant legislation. These quality management systems are constantly checked and recertified by certification bodies and public authorities. As a result, the company is constantly developing.

Internal, the quality is ensured by a quality management officer, a safety officer and a qualified person in accordance with MDR 2017/745.

The Funke Medical GmbH is distinguished by its outstanding medical competence, continuous high quality and reliable delivery capability.

Certified quality
management system
according to
DIN EN ISO 13485



Certified manufacturer
of CE-marked medical
devices according to
MDR 2017/745



Clinical evaluation according
to MEDDEV 2.7.1



Traceability of all products
with the help of the
UDI-Code



Laboratory studies:
microclimate, pressure relief,
flameretardancy,
biocompatibility



Quality management

qm@funke-medical.de

+49-2865-2669350



PRODUCT FEATURES





Clinic range

All Funke Medical decubitus therapy mattresses in the **clinic range** are characterized by their high quality and versatility.

Care bed mattresses
for all applications

- ✓ Normal ward
- ✓ Geriatric ward
- ✓ Intensive care unit
- ✓ Long-term care

High quality
decubitus therapy mattresses



Evacuation mattresses



- Your **all-in-one** solution -

More information on pages 29 – 32

BARIATRIC CARE



More information on pages 37 – 42

Climate-friendly mattresses



funke-medical.de/greenline

More information on pages 17 – 21

Hyper AIR® hybrid



Decubitus therapy air-cell mattress

More information on pages 43 and 46



Features of care bed mattresses

Clinic range



Efficient decubitus prophylaxis and therapy



High quality incontinence cover removable by zipper



CE-marking MDR 2017/745



Excellent Benefit



Easy cleaning



Flame retardant
DIN EN 597 (1+2)



High comfort for patients



Disinfection by wiping



Suitable for disinfection with chlorine *(max. 5000 ppm)*



Pain relief due to relieving positioning



Cover easy to disinfect



Biocompatible



Longer guarantee



Autoclavable



Maintenance-free



High stability and Long durability



Breathable open-cell material for an optimal microclimate



No risk of failure



Excellent pressure relief



Cover material virus- and bacteria-proof *accord. ASTM 1670 + 1671*



Economic



Ultrasonically welded seams



Waterproof but permeable to water vapour



Made in Germany



Firm edges



Bi-elastic cover to minimise shear forces



Climate-friendly



Zipper cover

Additional features of the **GREENLINE®** mattresses



Longer durability thanks to **ECO LAST®** foam



Longer guarantee



Approx. 22% less CO₂ per mattress

More information on
pages 17 – 21

Additional features of the **EVAQ®** mattresses



Special sliding material at the bottom



10 safe rescue handles for easier handling



Safety belts with quick-release buckle

More information on
pages 29 – 32



PRODUCTS



GREENLINE[©] — RESPONSIBLE CARE



— GREENLINE[®] —



RESPONSIBLE CARE

ECOLAST[®] Foam

Medical quality meets environment protection 4.0

approx. **22 %**
less CO₂

approx. 7,2 kg less CO₂ per mattress




Make your contribution
to CLIMATE PROTECTION!



GREENLINE[®] mattresses have a reduced CO₂ emission compared to conventional decubitus therapy mattresses (standard foam mattresses or alternating pressure mattresses).

*Source: Europur-Study 2015

Foam mattresses




1 kg foam = 3 kg CO₂

Conventional
foam mattresses




11 kg foam
= 33 kg CO₂

— GREENLINE[©] —




8,6 kg foam
= 25,8 kg CO₂

Alternating pressure mattresses



1 kg plastic = 3,5 kg CO₂



10 kg plastic
= 35 kg CO₂
+ rare earth elements
+ permanent energy consumption


ECOLAST® – Technical Advantages

GREENLINE® mattresses are made of ECOLAST® foam. Compared to conventional decubitus therapy mattresses (*standard foam mattresses or alternating pressure mattresses*) they have less CO₂ emission.


ECOLAST® foam is a specially developed hybrid foam. Compared to standard foam, the more open-cell structure ensures improved air circulation, so that the patient benefits from an optimal microclimate.

ECOLAST® foam offers perfect pressure relief and an improved durability.


Hardness loss after exposure to humidity
(35 °C, 95% RH*¹)



CO₂ emission per m³ in kg*³



Capability emission figure*⁴



*¹ RH = relative humidity

*² RG = density (kg/m³)

*³ Source: Europur-Study 2015

*⁴ In comparison:

Elasticity test ASTM D3574/H • Wet-Compression-Test ASTM D395 • Elongation ISO 1798 • Tensile strenght ISO 1798 • CO₂ emissivity

GREENLINE[©]




CLINIC CARE

— GREENLINE[®] —

RESPONSIBLE CARE

ECOLAST[®] Foam Technology



Hyper Foam PLUS *clinic*

More information about
GREENLINE[®] can be found
on pages 17 – 21

Your retailer:

PRS 
Pressure Relieving Solutions



Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

Applicable up to stage 3
(EPUAP)


Middle risk
(Norton-, Braden-, Waterlow-Scale)

Hyper Foam PLUS^{clinic}

— GREENLINE® —

Decubitus therapy mattress

5 YEARS GUARANTEE
*²



Size

Height 14 cm
Length See current price list
Width See current price list

GREENLINE®

ECOLAST® Foam technology:

- Ensures longer durability
- Less CO₂ emission



*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.


*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

— GREENLINE[®] —

RESPONSIBLE CARE

ECOLAST[®] Foam Technology



Hyper Foam 2 *clinic*

— GREENLINE[®] —

More information about
GREENLINE[®] can be found
on pages 17 – 21

Your retailer:

PRS 
Pressure Relieving Solutions



*approx. 7,2 kg less CO₂ per mattress (Source: Europur-Study 2015)

Funke Medical GmbH • Ährenfeld 10 • 46348 Raesfeld • Tel.: +49-2865-2669350 • Fax: +49-2865-8133 • info@funke-medical.de

klinikmatratze.de

-Germany-

Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

Applicable up to stage 4
(EPUAP)


Very high risk
(Norton-, Braden-, Waterlow-Scale)

Hyper Foam 2^{clinic}

GREENLINE®

Decubitus therapy mattress

7 YEARS GUARANTEE
*2



Size

Height 15 cm
Length See current price list
Width See current price list

GREENLINE®

ECOLAST® Foam technology:

- Ensures longer durability
- Less CO₂ emission



*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

Decubitus therapy mattresses with evacuation function



- Your **all-in-one** solution -





Decubitus therapy mattresses with evacuation function

- Your **all-in-one** solution-

In case of fire you have only **3 minutes** to escape!



All EVAQ® decubitus therapy systems are equipped with an evacuation function.
Patients can be evacuated by a single person in case of fire.

It takes less than **60 seconds** to evacuate one person from one floor to the other.

Watch our evacuation film and see for yourself.

Scan me! 



FILMLINK

Watch our evacuation film here!



- Your **all-in-one** solution -


Comparison – Furnishing of institutions

(e.g. clinics, elderly homes)

Furnishing without EVAQ®



+



Standard care mattress
mostly without
decubitus therapy function


+



Evacuation system
e.g. stretchers,
rescue sheets

In case of evacuation

+




4 – 6 rescuer

Furnishing with EVAQ®



+



evAQ®

=

mattress
+ decubitus therapy system
+ evacuation function

In case of evacuation

+



1 rescuer



Resource saving (money and workforce)!



Decubitus therapy mattresses with evacuation function

All mattresses in this range are equipped with an evacuation function. This means that you can evacuate patients quickly and easily without transferring them on rescue sheets or carry them on stretchers!


Evacuation systems such as stretchers and rescue sheets often require 2 or more persons to use them.

The EVAQ® decubitus therapy mattresses enable you to carry out evacuation gently using only one person, even via stairs!

§§ Quite apart from the fact that every facility operator is responsible for taking measures to ensure rapid evacuation! §§

SAFER → FASTER → GENTLER


- Safety belts with quick-release buckle
- Special sliding material at the bottom
- Stairway evacuation possible
- Firm edges provide stability
- Rescue handles on the cover



— GREENLINE[®] —

RESPONSIBLE CARE

ECOLAST[®] Foam Technology



evAQ[®] - PRO

— GREENLINE[®] —
- Your **all-in-one** solution -

More information about
GREENLINE[®] can be found
on pages 17 – 21

Your retailer:

PRS 
Pressure Relieving Solutions



Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning


Applicable up to stage 3
(EPUAP)

Middle risk
(Norton-, Braden-, Waterlow-Scale)

evAQ® - PRO
— GREENLINE® —

Decubitus therapy mattress
with evacuation function

5 YEARS GUARANTEE
*2



Size

Height 14 cm
Length See current price list
Width See current price list

GREENLINE®

ECOLAST® Foam technology:

- Ensures longer durability
- Less CO₂ emission



*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.


*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

— GREENLINE[®] —

RESPONSIBLE CARE

ECOLAST[®] Foam Technology



ev**AQ**[®] - PRO II

— GREENLINE[®] —

- Your **all-in-one** solution -

More information about
GREENLINE[®] can be found
on pages 17 – 21

Your retailer:

PRS 
Pressure Relieving Solutions



*approx. 7,2 kg less CO₂ per mattress (Source: Europur-Study 2015)

Funke Medical GmbH • Ährenfeld 10 • 46348 Raesfeld • Tel.: +49-2865-2669350 • Fax: +49-2865-8133 • info@funke-medical.de

klinikmatratze.de

-Germany-

Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

Applicable up to stage 4
(EPUAP)


Very high risk
(Norton-, Braden-, Waterlow-Scale)

evAQ® - PRO II

— GREENLINE® —

Decubitus therapy mattress
with evacuation function

7 YEARS GUARANTEE
*2



Size

Height 15 cm
Length See current price list
Width See current price list

GREENLINE®

ECOLAST® Foam technology:

- Ensures longer durability
- Less CO₂ emission



*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

GREENLINE[©]




BARIATRIC CARE

— GREENLINE[®] —

RESPONSIBLE CARE

ECOLAST[®] Foam Technology

Climate-friendly



MAXX 250
— GREENLINE[®] — *clinic*

More information about
GREENLINE[®] can be found
on pages 17 – 21

Your retailer:

PRS
Pressure Relieving Solutions



Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

Applicable up to stage 4
(EPUAP)


Very high risk
(Norton-, Braden-, Waterlow-Scale)

MAXX 250

— GREENLINE® — clinic

Bariatric decubitus therapy mattress

7 YEARS GUARANTEE *²



Size

Height 18 cm
Length See current price list
Width See current price list

GREENLINE®

ECOLAST® Foam technology:

- Ensures longer durability
- Less CO₂ emission



*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

— GREENLINE[®] —

RESPONSIBLE CARE

ECOLAST[®] Foam Technology



More information about
GREENLINE[®] can be found
on pages 17 – 21

More information about
EVAQ[®] can be found on
pages 29 – 32

evAQ[®] - PRO XL

— GREENLINE[®] —
- Your **all-in-one** solution -

Your retailer:

PRS
Pressure Relieving Solutions



Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

Applicable up to stage 4
(EPUAP)


Very high risk
(Norton-, Braden-, Waterlow-Scale)

evAQ[®] - PRO XL

— GREENLINE[®] —

Bariatric decubitus therapy mattress
with evacuation function

7 YEARS GUARANTEE
*2



Size

Height 18 cm
Length See current price list
Width See current price list

GREENLINE[®]

ECOLAST[®] Foam technology:

- Ensures longer durability
- Less CO₂ emission



*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!


*³ For more details see instructions for use.



AIRLINE[©]



SYSTEMATIC CLINIC CARE



Hyper AIR® hybrid - AIR-CELL MATTRESS -

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning


Applicable up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)

Hyper AIR[®] hybrid

Decubitus therapy air-cell mattress

5 YEARS
GUARANTEE
*²



Size

Height 16 cm
Length See current price list
Width See current price list

GREENLINE[®]

ECOLAST[®] Foam technology:

- Ensures longer durability
- Less CO₂ emission



Breathable
open-cell
Geltouch[®] foam +
ECOLAST[®] foam +
Visco foam



Autocla-
vable



Cover material
virus- and
bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Pressure
< 18 mmHg



Bi-elastic cover
to minimise
shear forces



CE- marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597
(I+II)




No risk of failure
due to self-inflating
FOAM-IN-CELL-
TECHNOLOGY



No electricity
required



Max. 5000 ppm



Patient weight:
0 – 250 kg

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.


*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

Hyper AIR[®] hybrid


- AIR-CELL MATTRESS -

WATCH NOW
on YouTube!



1

Put the mattress in the bed.
Please note the imprint on the top
of the mattress „THIS SIDE UP“.



2

Open the valve by pulling and wait
approx. 5 minutes, until the air-cells
have inflated automatically.

3

Close the valve
by pressing.

(Valve lock cannot fall off)

4

Place the patient on the mattress.
Carefully open the valve by pulling the
grey cap, so that air escapes and carry
out a fingertest* at the same time.

5

Close the valve by
pressing.


6

Daily checking of the correct filling
of the air-cells is recommended.
An adjustment may be necessary.


**Place hand between buttocks and mattress.
Optimal lying comfort is achieved when
only a slight finger movement of 1 – 2,5 cm
is possible.*



MEMOLINE[®]



SYSTEMATIC CLINIC CARE



Hyper Foam VISCO^{clinic}

Your retailer:

PRS
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning


Applicable up to stage 3
(EPUAP)

Middle risk
(Norton-, Braden-, Waterlow-Scale)

Hyper Foam VISCO^{clinic}

Decubitus therapy mattress

5 YEARS GUARANTEE
*²



Size

Height 14 cm
Length See current price list
Width See current price list




*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC CLINIC CARE



evAQ[®] - PRO VISCO

- Your **all-in-one** solution -

More information about
EVAQ[®] can be found on
pages 29 – 32

Your retailer:

PRS
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

Applicable up to stage 3
(EPUAP)

Middle risk
(Norton-, Braden-, Waterlow-Scale)

evAQ[®] - PRO VISCO

Decubitus therapy mattress
with evacuation function

5 YEARS
GUARANTEE

*2



Size

Height 14 cm
Length See current price list
Width See current price list



*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

POSITIONING AIDS



Positioning aids

POSIMED®

POSITIONING PILLOWS

-  Efficient decubitus prophylaxis and therapy
-  Pain relief due to relieving positioning
-  Rapid adaptation to the contours of the body



Icon image

POSIMED® HEEL PROTECT®

HEEL PROTECTOR BOOT

-  Efficient decubitus prophylaxis and therapy
-  Suitable for the treatment of diabetic foot ulcers




Icon image

MobiTurn®

TURNING AND POSITIONING AID

-  Enables back-friendly and effortless work
-  Easy patient turning



Icon image

POSIMED® Product Range

Universal Cushion



Can be used for all common positioning types

Size in cm	Item No.
37 / 26	PO-UNIK3726
56 / 40	PO-UNIG5640

Heel Rail



e.g. for positioning the heel

Size in cm	Item No.
39 / 21 / 8	PO-FERS3921

Forearm Splint



e.g. for positioning the forearm

Size in cm	Item No.
39 / 21 / 8	PO-UNTE3921

Abduction Pillow



e.g. for positioning the knee, ankle and epicondyle

Size in cm	Item No.
60 / 35	PO-SPRE6035

Bolster



e.g. for supporting and positioning the knee, neck and other parts of the body

Size in cm	Item No.
60 / 21	PO-ROLL6021

Hip Abduction Pillow



e.g. for positioning the hip – reduces the risk of dislocating hip prosthesis

Size in cm	Item No.
35 / 28 / 14	PO-HUEFT3528

Full Body Positioning Cushion



e.g. for support in the stable 30° position, with integrated pillow

Size in cm	Item No.
170 / 74	PO-LAGE0165
182 / 54	PO-LAGE1165

Semi-Circle



e.g. for upper body support

Size in cm	Item No.
137 / 80	PO-HALB1378

POSIMED®


Positioning pillows

POSIMED® pillows are positioning aids in many different shapes and sizes for pressure sore prevention and treatment support.


With the various positioning aids the patient can be moved marginally and the well-balanced pressure distribution enables the patient to lie in a relaxed and pressure-relieving position. The aim is to keep the duration and the intensity of pressure to a minimum. Our positioning system is based on a pressure-relieving and malleable ball-shaped filling.

Positioning examples


30° lateral position




Heel protection



Semi-sitting position



135° lateral position



SYSTEMATIC DECUBITUS THERAPY



Icon image

Your retailer:

POSIMED®
- POSITIONING CUSHION -
universal cushion

Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

For decubitus prophylaxis
and therapy

Cover is easy to disinfect

POSIMED®

- POSITIONING CUSHION -

universal cushion

5 YEARS
GUARANTEE
*2

Cushion:

- Washable up to 60°C
- Autoclavable
- Disinfection by wiping
- Disinfectant-proof*³
- Suitable for disinfection with chlorine*³

Biocompatible

Waterproof but breathable
incontinence cover ensures an
excellent microclimate

Cushion filling ensures:

- Rapid adaptation to the contours of the body
- Optimised pressure relief
- Physiological heat maintenance

To support the positioning of
extremities

Icon image

Item description

Size in cm

Item No.

Length x Width

Universal Cushion small

37 x 26

PO-UNIK3726

Universal Cushion large

56 x 40

PO-UNIG5640



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Optimised
pressure relief



Bi-elastic cover to
minimise shear
forces



CE marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)



Easy to clean by
wipe disinfection



Max. 5000 ppm

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties (pro rata basis)!

*³ For more details see instructions for use.

Static positioning cushions for the extremities help with decubitus prophylaxis and/or therapy. They are used to support pressure-free positioning and to change and stabilise the position of people being cared for, especially in bed. The positioning cushions cannot replace manual changes of position, but they may be able to extend the intervals between changes and facilitate positioning. The duration of the interval depends on the individual state of health of the person being cared for. To ensure the safe use of the positioning cushion, the carer must have the appropriate technical expertise or has been instructed in the use of the product by a care professional.

Indications:

- **Decubitus prophylaxis and therapy:** Gentle positioning allows vulnerable areas of the body such as arms, legs and/or heels to be temporarily positioned freely or softly.
- **Mobility restrictions and/or impaired body perception:** The positioning cushion helps to reduce the strength of pressure and shearing forces and to secure the position of people being cared for who are no longer able to change position independently.
- **Pain:** Pain relief through relieving positioning

Positioning examples

General preparation




A.

First check the universal cushion for cleanliness and possible damages.



B.

Please ensure that the filling is evenly distributed in the cushion during the positioning process.




C.

The care bed should be adjusted to the working height of the carer. The height of the carer's hip bone is used as a reference point.

If the person being cared for is responsive, they should be informed of every step. Make eye contact and, if accepted by the person being cared for, place your hand on their shoulder.


Positioning of the arm

The universal cushion is used for pressure-relieving positioning of vulnerable areas of the body such as arms.



1.


The person being cared for should be lying in a supine position. The carer is positioned at the arm that is to be relieved.



2.

The carer ensures not to lift the arm directly at the joints.

The carer lifts the forearm from underneath and with the other hand he places the cushion under the forearm.



3.

It is important that the elbow of the person being cared for is either exposed or resting completely on the cushion. The arm should not be lying on the edge.

If possible, make sure that the person being cared for agrees with the position and check the arm again to make sure it is well positioned.


To bring the person being cared for back into the initial position, the following measures should be taken:

The carer lifts the arm of the person to be cared for and removes the universal cushion. Check again that the person to be cared for is lying comfortably.

Positioning of the leg

The universal cushion is used for pressure-relieving positioning of vulnerable areas of the body such as legs and/or heels.

Prepare the positioning process as described in picture A-C.




The person being cared for should be lying in a supine position. The carer is positioned at the leg that is to be relieved.



The carer ensures not to touch the leg directly at the joints.

The carer lifts the patient's lower leg with one hand and places the universal cushion under the lower leg with the other hand.



The heel of the person being cared for should be exposed. If possible, make sure that the person being cared for agrees with the position and check the leg again for proper positioning.

To bring the person being cared for back into the initial position, the following measures should be taken:

The carer lifts the leg of the person to be cared for and removes the universal cushion. Check again that the person to be cared for is lying comfortably.

SYSTEMATIC DECUBITUS THERAPY



Icon image

Your retailer:

POSIMED®
- POSITIONING CUSHION- Bolster

Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

For decubitus prophylaxis
and therapy

Cover is easy to disinfect

POSIMED®

- POSITIONING CUSHION -

Bolster

5 YEARS
GUARANTEE
*2

Cushion:

- Washable up to 60°C
- Autoclavable
- Disinfection by wiping
- Disinfectant-proof*³
- Suitable for disinfection with chlorine*³

Biocompatible

Waterproof but breathable
incontinence cover ensures an
excellent microclimate

Cushion filling ensures:

- Rapid adaptation to the contours of the body
- Optimised pressure relief
- Physiological heat maintenance

Used for pressure-relieving
positioning of vulnerable areas
of the body such as the head,
neck and shoulder area, pelvis,
knees and legs.

Item description	Size in cm Length x Width	Item No.
Bolster	60 x 21	PO-ROLL6021



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Optimised
pressure relief



Bi-elastic cover to
minimise shear
forces



CE marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)



Easy to clean by
wipe disinfection



Max. 5000 ppm

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties (pro rata basis)!

*³ For more details see instructions for use.

Static positioning cushions for the partial body help with decubitus prophylaxis and/or therapy. They are used to support pressure-free positioning and to change and stabilise the position of people being cared for, especially in bed. The positioning cushions cannot replace manual changes of position, but they may be able to extend the intervals between changes and facilitate positioning. The duration of the interval depends on the individual state of health of the person being cared for. To ensure the safe use of the positioning cushion, the carer must have the appropriate technical expertise or has been instructed in the use of the product by a care professional.

Indications:

- **Decubitus prophylaxis and therapy:** Gentle positioning allows vulnerable areas of the body such as head, neck and shoulder area, pelvis, knees and legs to be temporarily positioned freely or softly.
- **Mobility restrictions and/or impaired body perception:** The positioning cushion helps to reduce the strength of pressure and shearing forces and to secure the position of people being cared for who are no longer able to change position independently.
- **Pain:** Pain relief through relieving positioning

Positioning examples

General preparation




A

First check the bolster cushion for cleanliness and possible damages.



B




C

The care bed should be adjusted to the working height of the carer. The height of the carer's hip bone is used as a reference point.

If the person being cared for is responsive, they should be informed of every step. Make eye contact and, if accepted by the person being cared for, place your hand on their shoulder.


Elevated positioning of the pelvis, knee and leg area

The bolster cushion is used for pressure-relieving positioning of vulnerable areas of the body such as the head, neck and shoulder area, pelvis, knees and legs.




1.

In this position, the bolster cushion is used to slightly bend the knee area. This reduces the tension on the abdominal wall and facilitates breathing. At the same time, a moderate amount of pressure is relieved in the sacrum/coccyx area. The pressure in this area is reduced as a hollow back is eliminated and the entire back serves as a support surface.



The carer places one arm under both lower legs of the person being cared for in order to carefully lift them. The carer can now use the other arm to place the bolster cushion underneath the lower part of the thighs.



The carer should ensure that the bolster cushion is not positioned directly in the hollow of the knee of the person being cared for in order to protect the nerves and blood vessels.

If possible, make sure that the person being cared for agrees with the position and check to ensure that they are positioned correctly.


To bring the person being cared for back into the initial position, the following measures should be taken:

The carer lifts the legs of the person to be cared for and removes the bolster cushion. Check again that the person to be cared for is lying comfortably.


Positioning of the head, neck and shoulder area

The bolster cushion is used for pressure-relieving positioning of vulnerable areas of the body such as the head.

Prepare the positioning process as described in picture A-C.



The person being cared for is positioned in the centre of the bed in a supine position.



With one hand, the carer lifts the cervical spine of the person being cared for.



With the other hand, the carer places the bolster cushion underneath the neck.



The carer should ensure that the bolster is not placed directly under the back of the head of the person being cared for.

If possible, make sure that the person being cared for agrees with the position and check the head again to make sure it is positioned correctly.

To bring the person being cared for back into the initial position, the following measures should be taken:

The carer lifts the cervical spine area of the person to be cared for and removes the bolster cushion. Check again that the person to be cared for is lying comfortably.

SYSTEMATIC DECUBITUS THERAPY



Icon Image

Your retailer:

POSIMED®

- POSITIONING CUSHION -

Full Body
Positioning Cushion

Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

For decubitus prophylaxis
and therapy

Cover is easy to disinfect

POSIMED®

- POSITIONING CUSHION -

Full Body
Positioning Cushion

5 YEARS
GUARANTEE
*²

Cushion:

- Washable up to 60°C
- Autoclavable
- Disinfection by wiping
- Disinfectant-proof*³
- Suitable for disinfection with chlorine*³

Biocompatible

Waterproof but breathable
incontinence cover ensures an
excellent microclimate

Cushion filling ensures:

- Rapid adaptation to the contours of the body
- Optimised pressure relief
- Physiological heat maintenance

For support in the
30° lateral and 135° lateral/
prone position

Item description

Size in cm

Item No.

Length x Width

Full Body Positioning Cushion Size M

170/74

PO-LAGE0165

Full Body Positioning Cushion Size L

182/54

PO-LAGE01165



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Optimised
pressure relief



Bi-elastic cover to
minimise shear
forces



CE marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)



Easy to clean by
wipe disinfection



Max. 5000 ppm

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties (pro rata basis)!

*³ For more details see instructions for use.

Static positioning cushions for the full body help with decubitus prophylaxis and/or therapy. They are used to support pressure-free positioning and to change and stabilise the position of people being cared for, especially in bed. The positioning cushions cannot replace manual changes of position, but they may be able to extend the intervals between changes and facilitate positioning. The duration of the interval depends on the individual state of health of the person being cared for. To ensure the safe use of the positioning cushion, the carer must have the appropriate technical expertise or has been instructed in the use of the product by a care professional.

Indications:

- **Decubitus prophylaxis and therapy:** Gentle positioning allows vulnerable areas of the the back of the head, shoulder blade, pelvis and especially the sacrum and coccyx.
- **Mobility restrictions and/or impaired body perception:** The positioning cushion helps to reduce the strength of pressure and shearing forces and to secure the position of people being cared for who are no longer able to change position independently.
- **Pain:** Pain relief through relieving positioning


Positioning examples

General preparation



A.

First check the full body positioning cushion for cleanliness and possible damages.



B.

Please ensure that the filling is evenly distributed in the cushion during the positioning process.




C.

The care bed should be adjusted to the working height of the carer. The height of the carer's hip bone is used as a reference point.

If the person being cared for is responsive, they should be informed of every step. Make eye contact and, if accepted by the person being cared for, place your hand on their shoulder.

Preparation for positioning

Bring the person to be cared for in the correct position so that the actual positioning can be carried out.



I.

If the person being cared for is lying in the centre or on the opposite side, they must be placed near the edge of the bed where the carer is standing.


To position the person to be cared for at the edge of the bed, the carer lifts the left shoulder with the right hand and slides the left hand under the shoulder blade up to the beginning of the spine.



II.


The carer's right hand then grabs the left shoulder of the person being cared for and turns it slightly so that it rests on the carer's left forearm.

When the left forearm is pulled out, the person being cared for turns onto their back and moves towards the carer. The person being cared for should now lie in front of the edge of the bed.



III.

This process is repeated with the hips according to the same principle. If the other side is positioned, this should be reversed.




IV.

Finally, the legs are raised slightly below the knees and laid down so that the person being cared for is lying straight again.

30° lateral position


The full body positioning cushion is used for pressure-relieving positioning of vulnerable areas of the body such as the back of the head, shoulder blade, pelvis and especially the sacrum and coccyx.

If the person to be cared for is incorrectly positioned in bed, please follow the instructions in Fig. I.–IV.



1.

The person being cared for should be in a supine position near the edge of the bed where the carer is standing. The carer stands on the side of the person to be cared for that needs to be relieved. The Full Body Positioning Cushion is placed on the edge of the bed next to the person being cared for, between the carer and the patient.



2.

The Full Body Positioning Cushion is positioned under one half of the body of the person being cared for. To do this, the carer lifts the left shoulder of the person to be cared for. The same procedure is now carried out with the hips.



3.

The legs are also being supported. With careful, controlled movements, the upper leg is placed on the Full Body Positioning Cushion to rest. Please ensure that you do not grasp the joints.



4.

If possible, make sure that the person being cared for agrees with the position and check the head and feet again to make sure they are positioned correctly.

To bring the person being cared for back into the initial position, the following measures should be taken:


The carer grabs the shoulder blade with the right hand and carefully lifts the person to be cared for. The upper part of the cushion can now be removed with the left hand. Proceed in the same way with the pelvic area. The Full Body Positioning Cushion can now be removed. The person to be cared for should then be repositioned in the centre of the bed. The procedure is the same as in steps I–IV. Check again that the person to be cared for is lying comfortably.

135° lateral/prone position

The full body positioning cushion is used for pressure-relieving positioning of vulnerable areas of the body such as the back of the head, shoulder blade, sacrum and coccyx.

Prepare the positioning process as described in picture A-C.

If the person to be cared for is incorrectly positioned in bed, please follow the instructions in Fig. I.–IV.




The carer stands on the side to be relieved of the person being cared for. The Full Body Positioning Cushion is placed on the opposite edge of the bed.



To secure the person being cared for, the side rails of the bed should be raised on the opposite side. The carer lifts the hip of the person to be cared for with the right hand and places the arm under the upper body with the palm of the hand facing upwards.




The outside leg of the person being cared for is positioned at an angle. The carer should take care not to grab the leg firmly, but only guide it with the right hand.



The carer now places one hand under the shoulder of the person to be cared for and the other hand under the sacrum. The person being cared for can now be turned towards the Full Body Positioning Cushion.

If more force is required, the knee of the person to be cared for can be used as an additional support.



If the person being cared for is in a prone position, the carer should check the position of the lower shoulder again and, if necessary, carefully pull the shoulder out backwards.



If possible, make sure that the person being cared for agrees with the position and check the head and feet again to make sure they are positioned correctly.

To bring the person being cared for back into the initial position, the following measures should be taken:

The carer grabs the shoulder of the person being cared for with one hand and the sacrum with the other hand. The person being cared for can be returned to the supine position by turning. Check again that the person to be cared for is lying comfortably.

SYSTEMATIC DECUBITUS THERAPY



Your retailer:

POSIMED®
- POSITIONIERUNGS CUSHION

Autoclavable according to
RKI*¹ guidelines

Pain relief due to
relieving positioning

For decubitus prophylaxis
and therapy

Cover is easy to disinfect

POSIMED®

Positionierungs Cushion

5 YEARS
GUARANTEE
*2

Cushion:

- Washable up to 60°C
- Autoclavable
- Disinfection by wiping
- Disinfectant-proof*³
- Suitable for disinfection with chlorine*³

Biocompatible

Waterproof but breathable
incontinence cover ensures an
excellent microclimate

Cushion filling ensures:

- Rapid adaptation to the contours of the body
- Optimised pressure relief
- Physiological heat maintenance

Item description	Size in cm	Item No.
	Length x Width	
Heel Rail	39 x 21	PO-FERS3921
Forearm Splint	39 x 21	PO-UNTE3921
Abduction Pillow	60 x 35	PO-SPRE6035
Hip Abduction Pillow	35 x 28	PO-HÜFT3528
Semi-Circle	137 x 80	PO-HALB1378

POSIMED® positioning cushions are available in many
different shapes and sizes.



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Optimised
pressure relief



Bi-elastic cover to
minimise shear
forces



CE marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)



Easy to clean by
wipe disinfection



Max. 5000 ppm

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties (pro rata basis)!


*³ For more details see instructions for use.

SYSTEMATIC DECUBITUS THERAPY



 **POSIMED®**
HEEL PROTECT®
- Heel protector boot -

Your retailer:



Autoclavable according to
RKI*¹ guidelines

Suitable for the treatment
of diabetic foot ulcers

For decubitus prophylaxis
and therapy

Cover is easy to disinfect

POSIMED® HEEL PROTECT®

Heel protector boot

1 YEAR
GUARANTEE*²

*²

Soft fastener

Soft fastener for
drop foot prophylaxis



Mesh³ fabric for
maximum comfort
and climate control

Icon image

Cover is easy to clean:

- Washable up to 95°C
- Suitable for tumble dryer
- Autoclavable
- Disinfection by wiping
- Disinfectant-proof*³
- Suitable for disinfection with chlorine*³

ROTATION-STOP



Unique
GELTOUCH®-technology
with cooling-effect

Size

Item No.

ONE SIZE FITS ALL

PO-HEELPR




Breathable
open-cell
Geltouch® foam +
HR foam



Autoclavable



Optimised
pressure relief



Bi-elastic cover
to minimise shear
forces



CE marking
MDR 2017/745



Designed to fit
a person of
any size




Max. 5000 ppm

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC POSITIONING



MobiTurn®

TURNING AND POSITIONING AID FOR
IMMOBILE PATIENTS

Your retailer:



Autoclavable according to
RKI*¹ guidelines

Enables back-friendly and
effortless work

Cover easy to disinfect

MobiTurn[®]


Turning and positioning aid
for immobile patients

2 YEARS
GUARANTEE
*2

Cover:

- Removable
- Washable up to 95°C
- Disinfection by wiping
- Disinfectant-proof*³
- Suitable for disinfection with chlorine*³
- Autoclavable

Ultrasonically welded
seams



Core: PU foam

Areas of application:

- Ergonomic shape ensures optimal patient positioning
- Structure enables safe and more painless turning of the patient
- Back-friendly and easy work for nurses and caring relatives



Cover material
virus- and bacte-
ria-proof accord.
ASTM 1670 + 1671



Waterproof but
permeable to
water vapour



Bi-elastic cover to
minimise shear
forces



CE-marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)



Autoclavable



Cover easy
to disinfect



Max. 5000 ppm



Back-friendly



Patient weight:
0- 250 kg



Easy patient
turning


*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties (pro rata basis)!

*³ For more details see instructions for use.



With the MobiTurn® positioning/turning aid, people being cared for can be turned safely and more painlessly from the supine to the side position in bed. The transfer by the caregiver is simplified and thus enables back-friendly and energy-saving work. In addition, the individual situation and potential illnesses must be taken into consideration when using the medical device. If you have any further questions, you should consider your doctor.



1.

The person being cared for lies flat in bed in the supine position. Please make sure that the person being cared for is not too close to the edge of the bed, if needed, the existing side rails should be raised on the side towards the person being cared for is to be turned. This prevents the person from falling out when being turned.



2.

The legs of the person being cared for are lifted by the care giver and the MobiTurn® is placed flat under the legs.




3.


The legs of the person being cared for need to be placed straight in the recesses of the MobiTurn® positioning/turning aid.



4.

In the next step, the MobiTurn® is straightened up and should be at the knees of the person to be cared for.






5.

Now the person to be cared for is prepared for turning. The arm facing away from the desired direction of positioning is placed on the abdomen of the person to be cared for. The other arm, which is on the side of the bed, is placed upside at a 90° angle.



6.


The carer places one hand on the MobiTurn® and the other hand on the shoulder of the person to be cared for in order to turn him into the side position as gently as possible.



7.

The patient is now in the lateral position and needs to be stabilized. Afterwards the MobiTurn® positioning/turning aid must be removed.

If the patient needs to be treated in the lateral position for a short time, the MobiTurn® is not removed until the patient has been returned to the supine position after the treatment.



8.

To return the patient to the supine position, the carer places one hand on the MobiTurn® at the hollow of the knee and the other hand on the shoulder of the patient and turns them back onto their back.

To remove the MobiTurn®, it is tilted towards the end of the bed until it lies flat on the mattress. The legs of the person being cared for are lifted and the MobiTurn® can easily be removed.



SERVICES





Services



EXPERT ADVICE

Certainly there are always questions that you would like to clarify with us personally on the phone or on site. For this purpose, a highly qualified advisory team from the fields of medicine, technology, finance and regulatory affairs is at your disposal.

phone: +49 (0)2865 266 93 50 | e-mail: info@funke-medical.de



FAQ

In our FAQ we answer the most frequently asked questions.

More information can be found on page 55.



FREE TRAINING

Qualified employees ensure successful decubitus prophylaxis and therapy. The Funke Medical Academy offers free, MDR-compliant certificate courses for decubitus therapy. Simply request your login data by e-mail.

More information can be found on page 53.



MAINTAINING THE QUALITY STANDARD

All our medical products meet the highest quality standards and are certified according to the MDR. This ensures better patient care and provides legal certainty while using our products.



Services

STRENGTHENING SUSTAINABILITY

Sustainability is not just an empty promise for us. In order to live up to our social and environmental responsibility, we are constantly working to make our production and products more sustainable.

More information can be found on page 8.

GUARANTEE OF THE ABILITY TO DELIVER

A secure supply chain with short delivery times plays a particularly important role in the medical sector. In order to be able to provide this, we produce exclusively in Germany with a capacity of 1,200 mattresses per day. In addition, we have a very large stock of finished goods to be able to absorb peaks in incoming orders.

REFERENCES

We would be pleased to inform you in a personal conversation about our many years of experience in supplying domestic and foreign hospitals with our high-quality decubitus therapy systems.





Free of charge!

Funke Medical Academy

Update your professional skills!



Do the course



Pass the test



Receive the
certificate



Done!

academy.funke-medical.de




HEALTH:ANGEL
WINNER
- 2019 -

„Digital further education in the health sector“



24 hours
Worldwide
All Devices



Funke Medical Academy

Update your professional skills!

According to national and international legislation, all professional groups in the health sector are obliged to be trained regularly.

The E-Academy of the Funke Medical GmbH offers this necessary further education as online courses with certificates.

- You can use the Funke Medical Academy to **train your own staff for free**.
- All courses are **up-to-date** with the latest scientific findings and are **constantly updated**.
- All courses have been developed by **pedagogues** with a scientific background.
- The Funke Medical E-Academy is and remains **free of charge**.



How does it work?

You receive a **registration key** and a **enrolment key** from the Funke Medical GmbH.

You and your employees can
register and log in with this keys at
academy.funke-medical.de

After completion of the training and successful execution of the online tests, you and your employees receive immediately a **certificate**, which can be used for submission and as verification document.



„We have had competition in the house for a long time,
why should we switch?“

WE ARE:

- ▶ Certified manufacturer of CE-marked medical devices according to MDR 2017/745
- ▶ Market leader in the field of decubitus prophylaxis and therapy for homecare in Germany

WE OFFER:

- ▶ Top quality at a reasonable price (no trading levels)
- ▶ Medical evidence
- ▶ Experiences in the international clinic market (Reference list)
- ▶ Own product development with highly qualified employees
- ▶ MDR-compliant products ensure legal certainty and traceability
- ▶ Secured supply chain: production in Germany, high storage capacity
- ▶ Very short delivery times
- ▶ Cost savings through intelligent production
- ▶ Improving sustainability through environmentally friendly products
- ▶ Competent medical and technical advice
- ▶ Free training opportunities

MEDICAL DIRECTOR / HYGIENE

What certificates does Funke Medical GmbH have?

- ▶ • CE-marking
- Clinical evaluation acc. to MEDDEV 2.7.1
- QM acc. to ISO 13485
- Laboratory studies: microclimate, pressure relief, biocompatibility, flame retardancy

How to clean the mattress core?

- ▶ By autoclaving (at least 50 minutes at 105 °C)

Does the foam suffer during autoclaving?

- ▶ No, we did a small test study to prove the opposite.

MEDICAL DIRECTOR / HYGIENE

How can the cover be cleaned?

- ▶ The cover can be washed at 95°C or cleaned by wiping disinfection according to the instructions for use.

Are the covers sewn?

- ▶ No, these are ultrasonically welded.
(Cover material virus and bacteria-proof acc. ASTM 1670+1671)



NURSING MANAGEMENT

Can the mattress really be used up to the advertised level of decubitus?

- ▶ Yes, because there are independent test procedures and clinical evaluations to confirm that.

Is visco foam the best?

- ▶ This is an advertising myth without clinical evidence. The fact is that viscose foam is closed-pored and has a negative effect to the microclimate. This then of course has a direct influence on the favouring of a decubitus. The claim that visco foam offers better pressure relief is also technically not measurable and therefore untrue. There are foams with significantly better measurement results.

How can a virus-proof cover be breathable?

- ▶ The PU laminate is a semi-permeable membrane.

Does a higher density increase durability?

- ▶ This is a myth that does not apply to modern raw materials in the 21st century. There are test methods that clearly disprove it.

Are the mattresses also suitable for heavy patients?

- ▶ Yes, the standard is 140 – 150 kg, but we also have mattresses with up to 250 kg patient weight (e.g.: MAXX 250 clinic GREENLINE).

PURCHASING DEPARTMENT

Is „Made in Germany / Europe“ too expensive?

- ▶ We will gladly provide proof of the opposite with the same quality!

Do the mattresses have a guarantee?

- ▶ Yes, depending on the model between 3 – 7 years on the functional properties of the mattress core.

Do you also offer free tests?

- ▶ Yes, of course.



PURCHASING / QUALITY MANAGEMENT

What references can Funke Medical GmbH provide?

- ▶ The reference list is available on request.

What certificates does Funke Medical GmbH?

- ▶ • CE-marking
- Clinical evaluation acc. to MEDDEV 2.7.1
- QM acc. to ISO 13485
- Laboratory studies: microclimate, pressure relief, biocompatibility, flame retardancy

What is your best-selling mattress?

- ▶ Hyper Foam 2- more than 1.000.000 pieces!

What are your capacities?

- ▶ We have a capacity of approx. 1,200 mattresses per day.

How long are the delivery times?

- ▶ For smaller quantities 48 hours for stock goods, for larger quantities 14 days after receipt of order.

BED CONTROL CENTER

For handling the mattresses, can they also be equipped with handles?

- ▶ Yes, but they are sewn on, which is technically impossible otherwise.

How often are the mattress cores cleanable / autoclavable?

- ▶ At least 100 times

How often can the covers be washed according to the instructions in the user manual?

- ▶ At least 100 times

If you have any further questions about our products or our service, please do not hesitate to contact us:

+49 (0)2865 266 93 50 | info@funke-medical.de



klinikmatratze.de

Via our Funke Medical social media channels
don't miss any more of our news:

